

PHASE 02

Sputtering-Annealing

QUICK OPERATOR GUIDE
Release 0.9.36 · Build 2026.08.11-r15

PURPOSE

Run the automated COSCON sputtering and oven-annealing cycles while the operator performs the physical argon leak-valve actions requested by the GUI.

WORKFLOW AT A GLANCE

1 · PREPARE Degas, or verified continuation	2 · ARGON Open valve only when prompted	3 · SPUTTER COSCON Operate and output verification	4 · SAFE STATE Standby, then close valve	5 · ANNEAL Oven hold; repeat configured cycles
--	--	---	---	---

BEFORE STARTING ✓ Keep other COSCON control software closed. ✓ Confirm pressure telemetry, oven PID communication and leak-valve position. ✓ Use Start without initial Degas only for a legitimate continuation run.	KEY FACTS <table><tr><td>Default cycles</td><td>3</td></tr><tr><td>Sputtering</td><td>20 min per cycle</td></tr><tr><td>Annealing</td><td>620 C for 10 min</td></tr><tr><td>Pressure limits</td><td>Warning 5e-5 mbar; emergency 1e-4 mbar</td></tr></table>	Default cycles	3	Sputtering	20 min per cycle	Annealing	620 C for 10 min	Pressure limits	Warning 5e-5 mbar; emergency 1e-4 mbar
Default cycles	3								
Sputtering	20 min per cycle								
Annealing	620 C for 10 min								
Pressure limits	Warning 5e-5 mbar; emergency 1e-4 mbar								
CONTROLS AND END STATES Prompted action Only perform the leak-valve action currently shown by the GUI. Abort / Safe Stop Requests COSCON safe state and oven PID 0 C, then closes Phase 02 only.	RUN DATA Saved under: Data Samples/Sputtering-Annealing Data/								

SAFETY - VERIFY THE REAL HARDWARE

COSCON Operate requires valid pressure, interlock and measured output. One specific transient Energy Overload may receive one controlled retry; a repeated fault aborts activation.